

Light. Measurement. Excellence.



TLS120Xe

High Power
Tuneable Light Source
(280-1100nm)



A high power monochromatic probe is the key to access the optical properties of components, materials and photodetectors.

- ▶ Fluorescence
- ▶ Photoluminescence
- ▶ Photochemistry
- ▶ Detector calibration
- ▶ Image sensor QE
- ▶ Deposition monitoring
- ▶ Fluorescence microscopy
- ▶ Reflectance
- ▶ Transmittance
- ▶ Absorbance

The TLS120Xe is a high power tuneable xenon light source that puts monochromatic light at your fingertips.



Monochromatic Light At Your Fingertips
Choose the wavelength you want, with direct front panel control or automation over USB.


Integrated Solution
Housed in a 19" unit, the TLS120Xe is easy to deploy and ideally suited to OEM applications.

High Power
A long-life, ultra-quiet 75W short-arc xenon light source ensures high optical power with excellent stability.


Continuous Tuning
A concave grating monochromator provides narrow bandwidth and high wavelength agility, giving you the freedom to select the wavelength you require.

The TLS120Xe is customised to your requirements with a choice of lamps and output coupling. Compatible with our portfolio of measurement accessories, including integrating spheres, collimators and relay optics, let the TLS120Xe deliver light your way.


1. Optical power at 5nm bandwidth with 1mm core SMA patchcord



2. Optical power at 20nm bandwidth with 3mm liquid light guide



Designed around a single benchtop enclosure, the TLS120Xe requires only a mains power input for operation. Adapting to a wide range of applications in research, industry and OEM, benefit from the power of monochromatic light at your fingertips.



1

Integrated Solution

Housing lamp, power supply and monochromator in a 19" benchtop enclosure, convenient to use and easy to deploy.

2

Wavelength Selection

Enjoy plug and play functionality, with direct front panel control or automation over USB.

3

Output Optics

Ensures efficient coupling to the fibre of choice and the option of beam attenuation using an iris iris diaphragm.

4

Display

An OLED display reports wavelength, bandwidth, burn time and stability.

5

Monitor Port

Ensure optical alignment of replacement lamps using the monitor port to record transmitted optical power.

6

Automation

Take full control of the TLS120Xe over USB using SCPI commands or using the spectral acquisition application, BenWin+.

7


Monochromatic Light Source

An ultra-quiet 75W short-arc xenon lamp feeds a 120mm focal length monochromator, ensuring high optical power and continuous wavelength tuning.

8

Coupling Options

A range of fibres and light guides are on offer to adapt to the spectral range in your application.



Light Source

Lamp type	Short-arc OFR xenon lamp	
Operating current	Constant current	
Coupling	Ellipsoidal reflector, 60mm diameter	
Lamp option	75W lamp	100W lamp
Operating current	5.4A	7.2A
Nominal lamp power and voltage	75W, 15V	100W, 14V
Nominal life time	2000 hours	500 hours

Opto-Mechanical

Configuration	Constant Deviation Angle
Focal length	120mm effective
Grating mount	Single grating on-axis turret
Grating type	Concave blazed holographic
Grating line density	1200
Nominal blaze wavelength	380nm
Slit type	Fixed Slit
Drive type	Stepper motor
Drive resolution	0.5 arcsec per step
Maximum drive speed	115° per sec

Optical Performance

Wavelength range	280nm – 1100nm
Wavelength step	0.0035nm
Bandwidth at FWHM with 0.74mm slit	5nm
Bandwidth at FWHM with 1.48mm slit	10nm
Bandwidth at FWHM with 2.96mm slit	20nm
Bandwidth at FWHM with 5.92mm slit	40nm

General

Interface	USB 2.0
Software control	BenWin+ spectral acquisition application, SCPI
Overall dimensions	300L x 460W x 185H (mm)
Weight	8.8 kg
Power supply	100-240V AC 50-60 Hz
Orientation	Horizontal only
Operating temperature range	10-35°C
Operating humidity range	30 % to 70 % (no condensation, less than 70 % above 30 °C)

Part Number	Description
TLS120Xe	High Power Tuneable Light Source (280nm-1100m)
19719	100W short-arc xenon lamp, nominal 500 hours lifetime
19583	75W short-arc xenon lamp, nominal 2000 hours lifetime
2x SA-FSP (0.72)	Fixed slit plate pair to set 5nm nominal bandwidth
2x SA-FSP (1.45)	Fixed slit plate pair to set 10nm nominal bandwidth
2x SA-FSP (2.90)	Fixed slit plate pair to set 20nm nominal bandwidth
2x SA-FSP (5.80)	Fixed slit plate pair to set 40nm nominal bandwidth
SMA-1500-1000	High OH ⁻ SMA fibre (250-1200nm), 1m long, 1.5mm diameter
FOP-UV-4-4-1000	UV grade fused silica fibre bundle (200-1350nm) 1 meter long, 4mm diameter
LLG-3-1000	Liquid light guide (280-700nm), 1m long, 3mm diameter

UK Contact Information

Bentham Instruments Limited
2 Boulton Road
Reading
Berkshire
RG2 0NH
United Kingdom

+44 118 975 1355
sales@bentham.co.uk
www.bentham.co.uk



Distributor Contact Information

AUSTRALIA

Photometric Solutions Intl Pty Ltd
sales@photometricsolutions.com

AUSTRIA

LOT-Quantum Design GmbH
info@lot-qd.de

BELARUS

UE PROFCON
marketing@profcon.by

CHINA

IVT Solar Pte Ltd
sales@ivt-solar.com
Sensing Instruments Co., Ltd,
sales@sensingm.com

FINLAND

Oy Mitaten Finland Ab
sales@mitaten.fi

FRANCE

Trioptics GmbH
tristan.athanaze@trioptics.fr

GERMANY

LOT-Quantum Design GmbH
info@lot-qd.de

INDIA

ATOS Instruments Marketing Services
atos@atosindia.com

INDONESIA

Industrial Vision Technology (S) Pte Ltd
wangyc@visiontec.com.sg

ISRAEL

IL Photonics
moshe@ilphotonics.com

ITALY

2M strumenti
info@2mstrumenti.com

JAPAN

Soma Optics Ltd
sales@somaopt.co.jp

KOREA

Won Woo Systems Co. Ltd
sales@wonwoosystem.co.kr

MALAYSIA

Industrial Vision Technology (S) Pte Ltd
wangyc@visiontec.com.sg

NEW ZEALAND

Photometric Solutions Intl Pty Ltd
sales@photometricsolutions.com

PORTUGAL

Lasing SA
info@lasing.com

RUSSIA

Ultratherm LLC
info@ultratherm.ru

SINGAPORE

Industrial Vision Technology (S) Pte Ltd
wangyc@visiontec.com.sg

SPAIN

Lasing SA
info@lasing.com

SWITZERLAND

LOT-Quantum Design AG
dumouchel@lot-qd.ch

TAIWAN

Asia Qtech Instrument Inc
service@qtechinstrument.com

U.S.A.

Market Tech, Inc.
info@markettechinc.net